

Nike Premier Lead Glasses



PREMIER-05
Matte Blue



PREMIER-07
Matte Black



PREMIER-08
Matte Anthracite

Lens Width:
60mm



Bridge Width:
13mm



Frame Width: 145mm

Temple Length:
135mm



Lens Height:
39.3mm



Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit.

FRAME SIZE:	Universal
WEIGHT:	76 G
PRESCRIPTION:	Plano, Single Vision
COLORS:	Matte Black, Matte Anthracite
PROTECTION	.75mm LE
INCLUDES	Anti-fog spray, safety strap add-on, protective case