



Nike Traverse Lead Glasses



TRAVERSE-01
Matte Black - Volt



TRAVERSE-02
Matte Anthracite - Red



Bridge Width: 13mm



Frame Width: 133mm



Temple Length: 135mm



Lens Width: 70mm



Lens Height: 42mm

The Nike Traverse lead glasses provide a secure fit and sharp vision. Meets ANSI requirement for side shielding

Frame Size	Universal
Prescription	Plano
Colors	Matte Black, Matte Anthracite
Weight	78.4G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case