

# Premium Wraps Lead Glasses



**683213**  
Tortoise Shell

- Only radiation protection frames on the market with lead glass lenses on the side shields for superior protection and excellent front and peripheral vision.
- Spring hinges, formable side temple bars and a fully adjustable nose piece
- Large front lens area is ideal for regular prescriptions.
- Prescription compatible (for most Rx)
- Lightweight and comfortable
- Large lens area for maximum protection

Lens Width: 57mm

Bridge Width: 19mm

Temple:  
140mm



Lens Height:  
59mm

Frame Width:  
130mm



<b>Frame Size</b>	Universal
<b>Prescription</b>	Single Vision (-SV) Plano
<b>Colors</b>	683213 Tortoise Shell
<b>Weight</b>	90G
<b>Protection</b>	Front: .75 mm LE Side: .75 mm LE
<b>Includes</b>	Anti-fog spray, safety strap add-on, protective case