

Premium Wraps Lead Glasses Single Vision, Plano



683213
Tortoise Shell

- Only radiation protection frames on the market with lead glass lenses on the side shields for superior protection and excellent front and peripheral vision.
- Spring hinges, formable side temple bars and a fully adjustable nose piece
- Large front lens area is ideal for regular prescriptions.
- Prescription compatible (for most Rx)
- Lightweight and comfortable
- Large lens area for maximum protection

Lens Width: 57mm

Bridge Width: 19mm

Temple:
140mm



Lens Height:
59mm

Frame Width:
130mm



Frame Size	Universal
Prescription	Single Vision (-SV) Plano
Colors	Tortoise Shell
Weight	97.3 G
Protection	Front: .75 mm LE Side: .75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case