

Xerxes Lead Glasses



80081 Shiny Black

Lens Width: 59mm

Bridge Width: 18mm

Temple Length:
140mm



Effective Diameter:
59mm



Frame Width:
129mm



- Wrap-around shape that delivers more protection from radiation with nose pads for added comfort.
- Lightweight and comfortable
- Large lens area for maximum protection
- Non Rx

Frame Size	Universal
Prescription	Plano
Colors	80081 Shiny Black
Weight	71G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case