

Maxi-Flex 5 Finger Lead Gloves



683300-112 Gray



683300-108 Blue

The perfect way to insure your hands and lower forearms are fully protected by the most flexible, comfortable radiation protection gloves on the market.

Comfortable fit with maximum dexterity. Soft removable, washable and replaceable inner liners* reduce odors and permit frequent laundering. Expanded durable vinyl outer covers are removable for inspection of lead core and permit replacement when they become worn. The lead vinyl cores are precision molded for consistent 0.5 mm lead equivalency and long life.



Maxi-Flex glove liners feature genuine Cool Wear moisture-wicking material

- .50 mm LE
- Tough, durable vinyl outer coverings
- Soft, Flexible and Comfortable
- Replaceable outer covers

1 PAIR

LEADED VINYL

Model	Description
683306	5 Replacement Liners
683307-108	5 Finger Gloves Outer Covers - Blue
683307-112	5 Finger Gloves Outer Covers - Gray