

Revolution Radiation Reduction Gloves

REVOLUTION



Infab's Exclusive Revolution Radiation Protection gloves are designed to be the most comfortable, flexible, lightweight and protective gloves on the market.



Latex and Powder Free

Every pair of Revolution Radiation Protection Gloves is completely latex and powder free.

Tactile Grip

Revolution Radiation Protection Gloves come with tactile grip finger tips. The tactile grip helps with keeping a strong grasp on instruments during a long procedure.

Three Protection Levels

There are three protection levels for Infab's Revolution Radiation Protection gloves: Light, Medium, and Heavy. Choose what best fits your needs.

Direct Beam Attenuation

Type	Direct Beam Energy Level		
	60 kVp	80 kVp	100 kVp
LIGHT	40%	35%	30%
MEDIUM	68%	61%	54%
HEAVY	75%	69%	63%

Values can vary by + 5%. 2.5 mm Al - Total Filtering

- Latex-free
- Powder-free
- Tactile grip
- Sterile
- Maximum tactile sensitivity
- Anatomical fit for long use without stress

Restorable with EO Gas

Boxes of 5 available with each model.
Add "-B"

1 PAIR

POWDER FREE

LEADED RUBBER

STERILE

	Model	Description
Light Direct Beam Attenuation 60 kV 40% (±5%) 80 kV 35% (±5%) 100 kV 30% (±5%) 2.5 mm Al - Total Filtering	700165	Latex-Free Radiation Reduction Gloves/Size 6.5/Pair
	700170	Latex-Free Radiation Reduction Gloves/Size 7/Pair
	700175	Latex-Free Radiation Reduction Gloves/Size 7.5/Pair
	700180	Latex-Free Radiation Reduction Gloves/Size 8/Pair
	700185	Latex-Free Radiation Reduction Gloves/Size 8.5/Pair
	700190	Latex-Free Radiation Reduction Gloves/Size 9/Pair

	Model	Description
Medium Direct Beam Attenuation 60 kV 68% (±5%) 80 kV 61% (±5%) 100 kV 54% (±5%) 2.5 mm Al - Total Filtering	735065	Latex-Free Radiation Reduction Gloves/Size 6.5/Pair
	735070	Latex-Free Radiation Reduction Gloves/Size 7/Pair
	735075	Latex-Free Radiation Reduction Gloves/Size 7.5/Pair
	735080	Latex-Free Radiation Reduction Gloves/Size 8/Pair
	735085	Latex-Free Radiation Reduction Gloves/Size 8.5/Pair
	735090	Latex-Free Radiation Reduction Gloves/Size 9/Pair

	Model	Description
Heavy Direct Beam Attenuation 60 kV 75% (±5%) 80 kV 69% (±5%) 100 kV 63% (±5%) 2.5 mm Al - Total Filtering	775065	Latex-Free Radiation Reduction Gloves/Size 6.5/Pair
	775070	Latex-Free Radiation Reduction Gloves/Size 7/Pair
	775075	Latex-Free Radiation Reduction Gloves/Size 7.5/Pair
	775080	Latex-Free Radiation Reduction Gloves/Size 8/Pair
	775085	Latex-Free Radiation Reduction Gloves/Size 8.5/Pair
	775090	Latex-Free Radiation Reduction Gloves/Size 9/Pair