

# Revolution Maxi-Flex 5 Finger Lead Gloves

REVOLUTION



Shown: 683300-501 Red/Black



Maxi-Flex glove liners feature genuine Cool Wear moisture-wicking material

- .50mm LE
- Tough, durable vinyl outer coverings
- Soft, Flexible and Comfortable
- Removable, washable and replaceable Cool Wear moisture-wicking inner liners
- Replaceable outer covers

1 PAIR LEADED VINYL

The perfect way to insure your hands and lower forearms are fully protected by the most flexible, comfortable radiation protection gloves on the market.

Comfortable fit with maximum dexterity. Soft removable, washable and replaceable inner liners\* reduce odors and permit frequent laundering. Expanded durable vinyl outer covers are removable for inspection of lead core and permit replacement when they become worn. The lead vinyl cores are precision molded for consistent 0.5 mm lead equivalency and long life.

Revolution Maxi-Flex 5 Finger Lead Gloves come in exclusively in our Revolution fabric, the softest and most comfortable material available for radiation protection.

Model	Description
683300-700	5 Finger 15" Length/Pair - Black/Black Revolution Material
683300-501	5 Finger 15" Length/Pair - Red/Black Revolution Material
683300-502	5 Finger 15" Length/Pair - Fuchsia/Black Revolution Material
683300-503	5 Finger 15" Length/Pair - Taupe/Black Revolution Material
683300-504	5 Finger 15" Length/Pair - Pink/Black Revolution Material
683300-505	5 Finger 15" Length/Pair - Navy/Black Revolution Material
683300-506	5 Finger 15" Length/Pair - Teal/Black Revolution Material
683306	5 Finger Gloves Replacement Liners
683307-700	5 Finger Gloves Outer Covers - Black/Black Revolution Material
683307-501	5 Finger Gloves Outer Covers - Red/Black Revolution Material
683307-502	5 Finger Gloves Outer Covers - Fuchsia/Black Revolution Material
683307-503	5 Finger Gloves Outer Covers - Taupe/Black Revolution Material
683307-504	5 Finger Gloves Outer Covers - Pink/Black Revolution Material
683307-505	5 Finger Gloves Outer Covers - Navy/Black Revolution Material
683307-506	5 Finger Gloves Outer Covers - Teal/Black Revolution Material