

Complete Patient Protection Set



683406
Storage Rack



Available in any of Infab's standard or premium colors



683401
9" H x 11.5" W
Toddler Apron



683402
13" H x 13" W
Child's Apron



683403
14" H x 17.5" W
Average Adult Apron



683404
17" H x 21" W
Large Adult Apron

- INCLUDES:** Storage rack with 4 aprons (Toddler, Child, Average, Large)
- COLORS:** Classic Nylon Colors
- CORE:** Lightweight Lead, KIARMOR Bi-Layer Low-Lead, KIARMOR Bi-Layer Lead-Free
- PROTECTION:** .25mm LE, .35mm LE, .50mm LE
- SIZES:** Toddler, Child, Average Adult, Large Adult

- Set or order the individual aprons
- Four sizes, from infant to large adult
- Equivalent to .50mm Pb
- Infab's exclusive belt with buckle closure
- Space efficient rack
- Solid aluminum, powder coated rack
- Rack dimensions - 1" wide and 8" deep

Model	Description
683400	Storage Rack with 4 Aprons (Toddler, Child, Average Adult and Large Adult)
683401	9" H X 11.5" W Toddler Apron Only
683402	13" H X 13" W Child's Apron Only
683403	14" H X 17.5" W Average Adult Apron Only
683404	17" H X 21" W Large Adult Apron Only
683405	Storage Rack with 1 Average Adult Apron and 1 Large Adult Apron
683406	Storage Rack Only