


Thyroid Shield



Comfortable, easy to wear, inexpensive and sanitary single use disposable shields reduce up to 60% of the direct beam between 50 and 120 kVp. Half the price of Bismuth shields but every bit as effective at blocking radiation.

Barium sulfate is one of the least harmful and most effective radiation attenuators. This is the same material administered during GI procedures and is the sole attenuating ingredient in our new environmentally safe Smart Armor Radiation protection eye and thyroid shields. It is safe enough to dispose of in the regular trash with no special handling required.



- Minimum 60% Dose Reduction
- Box of 50

Model	Description
PP-TS	Smart Armor Thyroid Shield - Box of 50

 **BARIUM SULFATE SHIELDS ARE NON-LATEX AND 100% DISPOSABLE**

*The American Journal of Roentgenology (AJR, March 2011, Vol. 196:3 pp.611-615)