


# Magnetic Thyroid Collar

## MTC



- Comfortable fit
- Large coverage area
- Optional monogram
- Optional KTC attachment

**Dimensions** | Hook and Loop Closures & Magnetic Closures



Model Number	Description
<b>MTC</b>	Magnetic Thyroid Collar
<b>Size</b>	X-Small, Standard, Medium
<b>Colors</b>	All
<b>Core</b>	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
<b>Protection</b>	.25mm LE, .35 mm LE, .50mm LE
<b>Closure</b>	Magnet
<b>Fabric</b>	Revolution or Classic
<b>Monogram</b>	<b>Block</b> or <i>Script</i>

Model Number	Description
MTC-XS	Magnetic Thyroid Collar, Size X-Small, Magnetic Closure
MTC-S1	Magnetic Thyroid Collar, Size Standard, Magnetic Closure
MTC-M1	Magnetic Thyroid Collar, Size Medium, Magnetic Closure