

# Tilted Mobile Barrier

## 683461



- Powder coated steel frame
- Acrylic impact panels, enclosed in a precision TIG welded frame for years of trouble-free service, completely enclosing the leaded glass panel for increased durability and safety
- 2 high-quality hospital grade locking casters
- 1.50 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Ergonomical viewing area
- Shipping worldwide by ground, air or sea freight

Model Number	Description
<b>683461</b>	Tilted Mobile Barrier - Lead
<b>683461-AC</b>	Tilted Mobile Barrier - Acrylic
<b>Dimensions</b>	70" H x 34 1/4" W X 25" D
<b>Glass Size</b>	24" H x 30" W
<b>Glass Options</b>	Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE
<b>Frame Colors</b>	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
<b>Frame Options</b>	Powder Coated Steel Stainless Steel (Upgrade)
<b>Curtain Color</b>	Revolution Grey Custom colors available.
<b>Wheels</b>	4 swivel casters, 2 with brakes
<b>Crate Size</b>	73" x 30" x 43"
<b>Crate Weight</b>	355 lbs / 161.03 kg

